

OMNISCOP DREAM



STEPHANIX
MEDICAL IMAGING SOLUTIONS



Omniscop DReam, surgical C-arm with Flat Panel Detector, is suitable for both hip and aortic prosthesis procedures.

Any kind of surgery is possible thanks to our wide range of generators and flat detectors.



COMFORT OF USE

The ergonomic design of the C-arm offers appropriate opening and positioning, from bariatric to extremities surgeries. The total size of the system has been designed to optimize the space inside the operating theatre.



The motorised variable height of the C-arm, the manual movements and the lateral handles on both sides, ensure an optimal positioning of the system.

Each release lever is marked with a colour code corresponding to one movement.



OMNISCOP DREAM





FLAT PANEL DETECTOR
IGZO TFT FOR 12X12 INCHES
A-SI FOR 8X8 INCHES

Two available sizes to fit with the specificity of each surgery.



DOSE REDUCTION

Half-dose exam, programmable and easy to reach in one-click. Removable grid.



POSITIONING LASER

Located on both tube and detector sides, for more precise centring.




WIRELESS FOOTSWITCH (OPTION)

To replace the standard wired footswitch.



ADVANCED POST-PROCESSINGS

Adapted to each preprogrammed anatomical protocol. For a reproducibility of image quality.




DSA AND SPECIFIC MEASUREMENT TOOLS (OPTIONS)

To meet all interventional vascular examinations needs.

USE



RF, DX AND XA MODALITIES



LATEST GENERATION OF DIGITAL RADIOGRAPHY SOLUTIONS: INTEGRATED, DYNAMIC, PORTABLE, WIRED AND WIRELESS DETECTORS.

REMOTE CONTROLLED TABLES



X-RAY ROOMS



X-RAY MOBILE UNITS



SURGICAL C-ARMS



STEPHANIX
MEDICAL IMAGING SOLUTIONS

Head office : 10, rue Jean Moulin
Z.I. du Bayon - 42150 La Ricamarie - FRANCE
Tel. +33 4 77 47 81 60 - Fax: +33 4 77 37 55 19
export@stephanix.com

www.stephanix.com